

schülke -t



Neutralizer for the automated reprocessing of dental devices

thermodent® neutralizer

Our Plus

- based on citric acid
- prevents deposits and coverings
- high level of material compatibility
- phosphate-free

Application areas

Neutralizer for automated treatment of dental instruments after alkaline cleaning (thermodent® powder, thermodent® alka clean), prevents coverings and to optimize cleaning results. Suitable for all common washers/disinfectors (WD).

Instructions for use

Dosage: approx. 0,1 - 0,2 % (1 - 2 ml/l).
Dosing with machine integrated dosing pumps.



thermodent® neutralizer

Product data

Components:
citric acid.

Chemical-physical data

Color	colourless
Density	ca. 1,17 g/cm ³ / 20 °C
Flash point	Not applicable
Form	liquid
pH	ca. 2 / 20 °C / concentrate
Viscosity, dynamic	ca. 4 mPa*s / Method : ISO 3219

Information for order

Item	Delivery form	Item no.
thermodent neutralizer 5 l KA	1/Canister	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Accessories

Accessories	Item no.
Can key for 5 + 10 l	135810

Related Products

- thermodent® alka clean
- thermodent® clear
- thermodent® powder

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com. For individual questions: Customer Care
Phone: +49 40 52100-666
E-Mail: info@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.



A company of the Air Liquide-Group

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 (0) 40 - 52100 - 0
Fax +49 (0) 40 - 52100 - 318
www.schuelke.com
mail@schuelke.com

Schülke & Mayr UK Ltd.
1 Jenkin Road
GB-Sheffield S9 1 AT
Phone +44 (0) 1142 - 5435 - 00
Fax +44 (0) 1142 - 5435 - 01
www.schulke.co.uk
mail.uk@schuelke.com

Schülke & Mayr Ges.m.b.H.
Seidengasse 9
1070 Wien, Austria
Phone +43 (0) 1 - 5232501-0
Fax +43 (0) 1 - 5232501-60
www.schuelke.at
office.austria@schuelke.com